



# Cervix Biopsy Punch

## Technical Data Sheet



Measurements per pack of 25: weight 0.5kg. Dimensions: 25cm Width, 31cm Height, 18.5cm Depth.

### Materials:

Acetal resin

Note that no latex is present

### Blister Pack:

Supplied in a sterile blister pack.

Do not use if pack is damaged.

### Safety Data:

There are no components presenting significant hazards or requiring specific personal protection or first aid measures.

All clinical waste should be disposed of in accordance with Local Health Authority and Department of Health guidelines.

Active Life Storage:- 5 years after manufacture. There is no specific safety recommendation for handling. Storage should be in the original packaging in a cool dry place.

The product should only be used by trained medical staff for examination in a clean environment.

Single-use only.



0120



A UK manufacturer - committed to quality, service, trust, innovation and value

**Pelican Healthcare Ltd**

Cardiff Business Park, Cardiff. CF14 5WF

Freephone: 0800 013 1220 Facsimile: 029 2074 7001

e-mail: [contactus@pelicanhealthcare.co.uk](mailto:contactus@pelicanhealthcare.co.uk)

[www.pelicanhealthcare.co.uk](http://www.pelicanhealthcare.co.uk)

**pelican**  
HEALTHCARE  
an eakin company