

**“The seal means I have been able to lead a relatively normal life as it stopped my bag leaking as soon as I started using them.”**

eakin Cohesive® Seal user



**eakin®  
Cohesive Range**

**Absorbency matters.**

# eakin<sup>®</sup> Cohesive – Seals

## Features & Benefits



eakin Cohesive Slims<sup>®</sup>

### Your Everyday Seal

- For everyday skin protection and security
- For frequent pouch changes

**Preventer**



Small eakin Cohesive<sup>®</sup>

### The Extra Protector

- For maximum skin protection
- For those with high fluid content

**Preventer**

**Problem Solver**



Large eakin Cohesive<sup>®</sup>

### Greatest Coverage

- For covering a large area of peristoma skin
- For complex skin problems or uneven skin surfaces

**Problem Solver**



eakin Cohesive StomaWrap<sup>™</sup>

### Your Time Saver

- For large, oval or irregular shaped stomas
- For those with limited dexterity

**Preventer**

**Problem Solver**

## eakin® Cohesive Seals & Stoma Wrap



**Our highly absorptive seals prevent moisture-associated skin damage (MASD) by offering:**

1. Super absorption 4.0g/g
2. High protection against enzymes (Elastase & Lipase)
3. Combatting irritants plus easier to mould and stretch increasing pouch wear time

“The difference one small Cohesive® Seal has made to my life is unbelievable. It has given my confidence back.”

## eakin® Cohesive Slims



- For everyday protection and security
- For frequent pouch changes
- For mild skin irritation
- For use with pre-cut pouches to create exact fit
- Ideal for use with convex pouches

“I use the Cohesive SLIMS® and they are superb. No skin reaction at all and they really last.”

**839005** - 30 per pack

**eakin® Cohesive**  
Seal – Small



- To provide maximum protection
- For those with high fluid output
- To resolve leakage problems
- To extend pouch wear time

**eakin® Cohesive**  
Seal – Large



- Ideal for Urostomates
- To cover a large area of peristomal skin
- Helps prevent irritation from pouch flange
- For use with deep convexity

**eakin® Cohesive**  
Stoma Wrap



- Unique shape
- Quick and easy to fit
- Ideal for large or oval shaped stomas
- No need to warm, stretch or tear to correct size
- For those with limited dexterity

“I have had a stoma since 1954, my skin had been raw and bleeding ever since or until a few months ago when I started using Cohesive®. Now my skin is beginning to heal and I no longer leak.”

**839002** -  
Small - 30 per pack  
**839001** -  
Large - 10 per pack

“I found them very easy to use and appreciate the fact that I didn’t have to tear them at all. I saw a massive improvement with my sore skin even after just the first day of wearing them.”

**839006** - 10 per pack

## eakin® Cohesive Cohesive Paste



- Alcohol free
- Absorbs moisture, keeping it away from the skin
- Will not sting or dry out
- Reduces digestive enzyme activity, stopping them from breaking down the skin

“I love Cohesive® Paste, it has given me great comfort and it’s easy to use. Skin ulcers around my stoma have now healed.”

**839010** - 60g Tube

## eakin® Support Belt



- Designed to provide full and complete abdominal support for prevention and treatment of peristomal hernias
- Lightweight, breathable and washable fabric
- Neutral tone fabric and slim line fit

**Find codes overleaf**

\* Assessments and measurement for support belt should be carried out by a Stoma Care Nurse



# eakin® Cohesive

eakin Cohesive®			
Product	Pack Size	Size	Code
Small Cohesive® Skin Barriers	5	10x10cm	839004
Large Cohesive® Skin Barriers	5	10x20cm	839003
Cohesive® Paste	60g Tube	-	839010

eakin Cohesive®				
Product	Pack Size	Diameter	Thickness	Code
Large Cohesive® Seal	10	98mm	2.1mm	839001
Small Cohesive® Seal	30	48mm	4.2mm	839002
Cohesive SLIMS®	30	48mm	3.0mm	839005
Cohesive StomaWrap™	10	85mm	3.0mm	839006



eakin® Support Belts			
Description	Measurements	20cm Wide	26cm Wide
Small	75-90cm	839051	839061
Medium	90-105cm	839052	839062
Large	105-120cm	839053	839063
Extra Large	120-135cm	839054	839064