

“It stayed put even after riding my bike”



pelican  
platinum

Platinum

Experience  
comfort, security  
and flexibility.

“Pelican Platinum is far better than pouches I have tried before. I have confidence now, they are supple and stick better to my skin.”

Platinum Contour Closed User

# Pelican Platinum

## Features & Benefits

pelican  
platinum

### Skin Protector

- 6 Contoured Edges
- Foam backed
- Reduces leakages and increases security
- Flexibility and freedom
- Comfort



### Split Fabric Backing

- Ease of use
- Access to filter



### Coloured Measuring Guide & Thumb Tab

- Easy to see
- Improves cutting accuracy
- Easy to grasp



### Twist Tap & Bung Options

- Reduced risk of accidents



### Soft Convexity

- Comfort
- Flexibility



### Non-return Valve

- Reduced risk of infection
- Urostomy only


### Discreet Outlet

- Central position
- Velcro dot
- New easy open outlet



### Hydrophobic Filter

- Avoid filter leakages
- Reduced pancaking and ballooning



### Pebble Shape & Range of Sizes

- Pre-cut and cut to fit range of sizes
- 12-60mm in convex
- 12-80mm in flat



## Pelican Platinum

### Closed

Platinum Contour Closed			Pack of 30
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-65mm	207212	207412	207612
22.5mm	-	207422	-
25mm	-	207425	-
27.5mm	-	207427	-
30mm	-	207430	-
32.5mm	-	207432	-
35mm	-	207435	-
37.5mm	-	207437	-
40mm	-	207440	-
45mm	-	207445	-
50mm	-	207450	-
55mm	-	207455	-
Cut to fit 15-80mm	-	208415	-



Platinum Contour Closed



Platinum Closed

Platinum Closed		Pack of 30				
Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
Cut to fit 12-65mm	107212	107112	107412	107312	107612	107512
22.5mm	-	-	107422	107322	-	-
25mm	-	-	107425	107325	-	-
27.5mm	-	-	107427	107327	-	-
30mm	-	-	107430	107330	-	-
32.5mm	-	-	107432	107332	-	-
35mm	-	-	107435	107335	-	-
37.5mm	-	-	107437	107337	-	-
40mm	-	-	107440	107340	-	-
45mm	-	-	107445	107345	-	-
50mm	-	-	107450	107350	-	-
55mm	-	-	107455	107355	-	-
Cut to fit 15-80mm	-	-	108415	108315	-	-

## Pelican Platinum

### Convex Closed

Platinum Contour Convex Closed			Pack of 10
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	407212	407412	407612
20mm	-	407420	-
25mm	-	407425	-
30mm	-	407430	-
35mm	-	407435	-
40mm	-	407440	-
Cut to fit 15-60mm	-	408415	-



Platinum Contour Convex Closed

Platinum Convex Closed			Pack of 10
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	307212	307412	307612
20mm	-	307420	-
25mm	-	307425	-
30mm	-	307430	-
35mm	-	307435	-
40mm	-	307440	-
Cut to fit 15-60mm	-	308415	-



Platinum Convex Closed

pelican  
platinum

## Pelican Platinum Drainable

Platinum Contour Drainable			Pack of 30
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-65mm	218212	218412	218612
20mm	-	218420	-
22.5mm	-	218422	-
25mm	-	218425	-
27.5mm	-	218427	-
30mm	-	218430	-
32mm		218432	
35mm		218435	
37.5mm		218437	
40mm		218440	
45mm		218445	
50mm		218450	
55mm		218455	
Cut to fit 15-80mm	-	220415	-



Platinum Contour  
Drainable

## Pelican Platinum Convex Drainable

Platinum Contour Convex Drainable			Pack of 10
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	418212	418412	418612
20mm	-	418420	-
25mm	-	418425	-
30mm	-	418430	-
35mm	-	418435	-
40mm	-	418440	-
Cut to fit 15-60mm	-	420415	-



Platinum Contour  
Convex Drainable

Platinum Drainable			Pack of 30
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-65mm	118212	118412	118612
20mm	-	118420	-
22.5mm	-	118422	-
25mm	-	118425	-
27.5mm	-	118427	-
30mm	-	118430	-
32.5mm		118432	
35mm		118435	
37.5mm		118437	
40mm		118440	
45mm		118445	
50mm		118450	
55mm		118455	
Cut to fit 15-60mm	-	120415	-



Platinum Drainable

Platinum Convex Drainable			Pack of 10
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	318212	318412	318612
20mm	-	318420	-
25mm	-	318425	-
30mm	-	318430	-
35mm	-	318435	-
40mm	-	318440	-
Cut to fit 15-60mm	-	320415	-



Platinum Convex  
Drainable

pelican  
platinum

## Pelican Platinum Convex Urostomy

Platinum Contour Convex Urostomy (with Bung) Pack of 10			
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	421212	421412	421612
20mm	-	421420	-
25mm	-	421425	-
30mm	-	421430	-
35mm	-	421435	-
40mm	-	421440	-
Cut to fit 15-60mm	-	422415	-



Platinum Contour  
Convex Urostomy  
with Bung

Platinum Contour Convex Urostomy (with Tap) Pack of 10			
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	423212	423412	423612
20mm	-	423420	-
25mm	-	423425	-
30mm	-	423430	-
35mm	-	423435	-
40mm	-	423440	-
Cut to fit 15-60mm	-	424415	-



Platinum Contour  
Convex Urostomy  
with Tap

## Pelican Platinum Convex Urostomy

pelican  
platinum

Platinum Convex Urostomy (with Bung) Pack of 10			
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	321212	321412	321612
20mm	-	321420	-
25mm	-	321425	-
30mm	-	321430	-
35mm	-	321435	-
40mm	-	321440	-
Cut to fit 15-60mm	-	322415	-



Platinum Convex  
Urostomy

Platinum Convex Urostomy (with Tap) Pack of 10			
Sizes	Mini Clear	Std Clear	Maxi Clear
Cut to fit 12-40mm	323212	323412	323612
20mm	-	323420	-
25mm	-	323425	-
30mm	-	323430	-
35mm	-	323435	-
40mm	-	323440	-
Cut to fit 15-60mm	-	324415	-



Platinum Convex  
Urostomy with Tap