

"The baseplate is more flexible during my work as a gardener. I bend and stretch a lot and the product feels nice and does not fold."



pelican  
platinum  
2 Piece 

# 2 Piece

For when life just clicks

"Very secure all round.  
I felt very confident."

# Platinum 2 Piece Features & Benefits

pelican  
**platinum**  
2 Piece



## Purple Guide Ring

- Easier location for attachment

## Pop Up Back Plate

- For added finger space
- Quick and easy bag changes
- Reduced pressure on abdomen/stomach



## Unique Peel Down Action

- For greater control of bag changes
- Ability to burp the bag




## Audible Click

- Know bag is securely attached
- Bag attaches at any angle



## Discreet Outlet

- Central position
- Velcro dot
- New easy open outlet



## 6 Contoured Edges



## Coloured Measuring Guide & Thumb Tab

- Easy to see
- Improves cutting accuracy
- Easy to grasp

## Longer Weartime Hydrocolloid

- Foam backed for extra comfort

## Platinum 2 Piece 60mm

## Platinum 2 Piece 80mm



Platinum 2 Piece Base Plate		Pack of 5
Size	Code	Stoma Sizes
60mm	600015	Cut to fit 15-37.5mm
	600020	20mm
	600022	22.5mm
	600025	25mm
	600027	27.5mm
	600030	30mm
	600032	32.5mm
	600035	35mm
600037	37.5mm	



Platinum 2 Piece Base Plate		Pack of 5
Size	Code	Stoma Sizes
80mm	800015	Cut to fit 15-50mm
	800040	40mm
	800045	45mm
	800050	50mm



Platinum 2 Piece Closed Pouch		Pack of 30
Sizes	Code	Pouch Size
60mm	607215	Mini
	607415	Standard
	607615	Maxi



Platinum 2 Piece Closed Pouch		Pack of 30
Sizes	Code	Pouch Size
80mm	807415	Standard
	807615	Maxi



Platinum 2 Piece Drainable Pouch		Pack of 30
Sizes	Code	Pouch Size
60mm	618215	Mini
	618415	Standard
	618615	Maxi



Platinum 2 Piece Drainable Pouch		Pack of 30
Sizes	Code	Pouch Size
80mm	818415	Standard
	818615	Maxi

