

Product Guide

Ostomy Care Product Range



Contents

Style meets comfort



**Practical, comfortable, and stylish,
Pelican ModaVi ticks all the boxes:**

- ✓ Shape & Structure | ✓ Colour
- ✓ Fabric | ✓ Outlet | ✓ Filter
- ✓ Window | ✓ Fold

Pelican ModaVi



2PC Closed Bag	6
2PC Drainable Bag	6
Closed Bag	7
Soft Convex Closed Bag	8
Firm Convex Closed Bag	8
Drainable Bag	9
Soft Convex Drainable Bag	10
Firm Convex Drainable Bag	10

Features & benefits

Filter

Improved design can help to reduce pancaking and ballooning.

Shape & Structure

A new profile shape and weld, ensuring a more natural fit for a more comfortable experience.

Window

The large viewing window enables easy application & monitoring.

Colour

A choice of two stylish colours for more discretion.

Fold

Ability to change the size and shape of the bag, for a more discreet look.

Outlet

Improved outlet. Wider, simpler and easier to use when emptying the bag.

Fabric

New water-resistant fabric to keep patients feeling fresh for longer.



Pelican ModaVi 2PC Closed Bag

Units per box: 30

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
P2CS6N	P2CS6B	Pelican ModaVi 2 Piece Closed Bag	S	60mm Ring
P2CM6N	P2CM6B	Pelican ModaVi 2 Piece Closed Bag	M	60mm Ring
P2CL6N	P2CL6B	Pelican ModaVi 2 Piece Closed Bag	L	60mm Ring



Pelican ModaVi 2PC Drainable Bag

Units per box: 30

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
P2DS6N	P2DS6B	Pelican ModaVi 2 Piece Drainable Bag	S	60mm Ring
P2DM6N	P2DM6B	Pelican ModaVi 2 Piece Drainable Bag	M	60mm Ring
P2DL6N	P2DL6B	Pelican ModaVi 2 Piece Drainable Bag	L	60mm Ring



Pelican ModaVi Closed Bag

Units per box: 30

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
PCS65N	PCS65B	Pelican ModaVi Closed Bag	S	CTF 12-65mm
PCM65N	PCM65B	Pelican ModaVi Closed Bag	M	CTF 12-65mm
PCL65N	PCL65B	Pelican ModaVi Closed Bag	L	CTF 12-65mm
PCL80N	PCL80B	Pelican ModaVi Closed Bag	L	CTF 15-80mm
PCM25PN	PCM25PB	Pelican ModaVi Closed Bag	M	Pre-Cut 25mm
PCM30PN	PCM30PB	Pelican ModaVi Closed Bag	M	Pre-Cut 30mm
PCM35PN	PCM35PB	Pelican ModaVi Closed Bag	M	Pre-Cut 35mm
PCM40PN	PCM40PB	Pelican ModaVi Closed Bag	M	Pre-Cut 40mm



The Pelican Podcast

Let's Talk... dives into the lives of inspiring ostomates as they open up about living with a stoma. Hosted by Louise, an ostomate herself and a Pelican ambassador, who navigates a range of stoma-related topics with other ostomates, who share their stories with the aim of helping those who may feel isolated, alone or not understood during their stoma journey.



Listen to the podcast by scanning the QR code or search "Pelican Let's Talk..." on your chosen Podcast platform.

Pelican ModaVi Soft Convex Closed Bag

Units per box: 10

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
PCSS40N	PCSS40B	Pelican ModaVi Soft Convex Closed Bag	S	CTF 12-40mm
PCSM40N	PCSM40B	Pelican ModaVi Soft Convex Closed Bag	M	CTF 12-40mm
PCSL40N	PCSL40B	Pelican ModaVi Soft Convex Closed Bag	L	CTF 12-40mm
PCSL60N	PCSL60B	Pelican ModaVi Soft Convex Closed Bag	L	CTF 15-60mm
PCSM25PN	PCSM25PB	Pelican ModaVi Soft Convex Closed Bag	M	Pre-Cut 25mm
PCSM30PN	PCSM30PB	Pelican ModaVi Soft Convex Closed Bag	M	Pre-Cut 30mm
PCSM35PN	PCSM35PB	Pelican ModaVi Soft Convex Closed Bag	M	Pre-Cut 35mm
PCSM40PN	PCSM40PB	Pelican ModaVi Soft Convex Closed Bag	M	Pre-Cut 40mm



Pelican ModaVi Firm Convex Closed Bag

Units per box: 10

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
PCFS40N	PCFS40B	Pelican ModaVi Firm Convex Closed Bag	S	CTF 12-40mm
PCFM40N	PCFM40B	Pelican ModaVi Firm Convex Closed Bag	M	CTF 12-40mm
PCFL40N	PCFL40B	Pelican ModaVi Firm Convex Closed Bag	L	CTF 12-40mm
PCFL60N	PCFL60B	Pelican ModaVi Firm Convex Closed Bag	L	CTF 15-60mm



Pelican ModaVi Drainable Bag

Units per box: 30

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
PDS65N	PDS65B	Pelican ModaVi Drainable Bag	S	CTF 12-65mm
PDM65N	PDM65B	Pelican ModaVi Drainable Bag	M	CTF 12-65mm
PDL65N	PDL65B	Pelican ModaVi Drainable Bag	L	CTF 12-65mm
PDL80N	PDL80B	Pelican ModaVi Drainable Bag	L	CTF 15-80mm
PDM25PN	PDM25PB	Pelican ModaVi Drainable Bag	M	Pre-Cut 25mm
PDM30PN	PDM30PB	Pelican ModaVi Drainable Bag	M	Pre-Cut 30mm
PDM35PN	PDM35PB	Pelican ModaVi Drainable Bag	M	Pre-Cut 35mm
PDM40PN	PDM40PB	Pelican ModaVi Drainable Bag	M	Pre-Cut 40mm



Here's what our ambassadors had to say...



"After trialling Pelican ModaVi, the decision to change from my previous bag was a no-brainer. The bag looks good, feels good and performs outstandingly. Whether at work, rest or play, ModaVi won't let you down."

Keith



"The design is amazing, the shape and features are exactly what I want in an ostomy bag. The feel of the bag is something else, it's so soft against my skin."

Ant



"What magic have you put on this bag? I would 100% recommend it to a friend, to anyone, to the world!"

Rachel

Pelican ModaVi Soft Convex Drainable Bag

Units per box: 10

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
PDSS40N	PDSS40B	Pelican ModaVi Soft Convex Drainable Bag	S	CTF 12-40mm
PDSM40N	PDSM40B	Pelican ModaVi Soft Convex Drainable Bag	M	CTF 12-40mm
PDSL40N	PDSL40B	Pelican ModaVi Soft Convex Drainable Bag	L	CTF 12-40mm
PDSL60N	PDSL60B	Pelican ModaVi Soft Convex Drainable Bag	L	CTF 15-60mm
PDSM25PN	PDSM25PB	Pelican ModaVi Soft Convex Drainable Bag	M	Pre-Cut 25mm
PDSM30PN	PDSM30PB	Pelican ModaVi Soft Convex Drainable Bag	M	Pre-Cut 30mm
PDSM35PN	PDSM35PB	Pelican ModaVi Soft Convex Drainable Bag	M	Pre-Cut 35mm
PDSM40PN	PDSM40PB	Pelican ModaVi Soft Convex Drainable Bag	M	Pre-Cut 40mm



Pelican ModaVi Firm Convex Drainable Bag

Units per box: 10

Neutral Fabric Product Code	Black Fabric Product Code	Product Description	Size	Type
PDFS40N	PDFS40B	Pelican ModaVi Firm Convex Drainable Bag	S	CTF 12-40mm
PDFM40N	PDFM40B	Pelican ModaVi Firm Convex Drainable Bag	M	CTF 12-40mm
PDFL40N	PDFL40B	Pelican ModaVi Firm Convex Drainable Bag	L	CTF 12-40mm
PDFL60N	PDFL60B	Pelican ModaVi Firm Convex Drainable Bag	L	CTF 15-60mm





Choose freedom, choose eakin freeseal[®]

Freedom from leaks and sore skin
to live life with confidence.

Key Features



1.8mm thin

Our thinnest seal allows for minimal impact on the fit of the bag. Its thinness helps keep the bag discreet, whilst giving peace of mind against leaks.



Absorbent

eakin freeseal[®] is an absorbent seal, making it effective at stopping leaks.¹ eakin freeseal[®] actively absorbs output preventing leaks and therefore protects against skin damage.



Everyday Essential

eakin freeseal[®] is easy to use. It is flexible and mouldable allowing it to fit securely and comfortably to the shape of the stoma. It is ideal for use as part of a regular stoma care routine.



Easy to remove

eakin freeseal[®] is easy to remove, retaining its form to peel away in one piece. This allows for quicker and simpler bag changes.

"I've been using it since it came out and it's revolutionised my stoma bag. It's so easy to shape to fit around my stoma and once it's on I can't even feel it."
— **Summer, eakin freeseal[®] user**



Scan me
to request
a sample

1. Mc Groggan G, Haughey S and McDowell K (2018) An absorbent, enzyme-inhibiting seal reduces peristomal skin complications. *Gastrointestinal Nursing* 16:1 42-4

Platinum 2-Piece



60mm Baseplates	14
Soft Convex 60mm Baseplates	14
60mm Bags	15
80mm Baseplates	16
Soft Convex 80mm Baseplates	16
80mm Bags	16

Features & benefits

Purple guide ring

A visual aid for easier location for attachment. Also allows the bag to be attached at any angle.

Pop up back plate

For quick and easy bag changes that reduces pressure on the abdomen.

6 contoured edges

Adapts to any body shape for flexibility and comfort — particularly useful for hernias.

Unique peel down action

For greater control of bag changes.

Audible click

Gives confidence that the bag is secure.



Platinum 2PC 60mm Ring Size Baseplates

Units per box: 5

Product Code	Product Description	Size	Type
600015	Platinum 2 Piece Flat Baseplate	CTF 15-37.5mm	60mm Ring
600020	Platinum 2 Piece Flat Baseplate	Pre-Cut 20mm	60mm Ring
600022	Platinum 2 Piece Flat Baseplate	Pre-Cut 22.5mm	60mm Ring
600025	Platinum 2 Piece Flat Baseplate	Pre-Cut 25mm	60mm Ring
600030	Platinum 2 Piece Flat Baseplate	Pre-Cut 30mm	60mm Ring
600032	Platinum 2 Piece Flat Baseplate	Pre-Cut 32.5mm	60mm Ring
600035	Platinum 2 Piece Flat Baseplate	Pre-Cut 35mm	60mm Ring
600037	Platinum 2 Piece Flat Baseplate	Pre-Cut 37.5mm	60mm Ring



Platinum 2PC Soft Convex 60mm Ring Size Baseplates

Units per box: 30

Product Code	Product Description	Size	Type
600015C	Platinum 2 Piece Soft Convex Baseplate	CTF 15-40mm	60mm Ring
600020C	Platinum 2 Piece Soft Convex Baseplate	Pre-Cut 20mm	60mm Ring
600025C	Platinum 2 Piece Soft Convex Baseplate	Pre-Cut 25mm	60mm Ring
600030C	Platinum 2 Piece Soft Convex Baseplate	Pre-Cut 30mm	60mm Ring
600035C	Platinum 2 Piece Soft Convex Baseplate	Pre-Cut 35mm	60mm Ring
600040C	Platinum 2 Piece Soft Convex Baseplate	Pre-Cut 40mm	60mm Ring



Platinum 2PC 60mm Ring Size Closed Bags

Units per box: 30


Product Code	Product Description	Size
607215	Platinum 2 Piece Closed Bag	Mini
607415	Platinum 2 Piece Closed Bag	Standard
607615	Platinum 2 Piece Closed Bag	Maxi



Platinum 2PC 60mm Ring Size Drainable Bags

Units per box: 30


Product Code	Product Description	Size
618215	Platinum 2 Piece Drainable Bag	Mini
618415	Platinum 2 Piece Drainable Bag	Standard
618615	Platinum 2 Piece Drainable Bag	Maxi



Platinum 2PC 80mm Ring Size Baseplates

Units per box: 5


Product Code	Product Description	Size	Type
800015	Platinum 2 Piece Flat Baseplate	CTF 15-50mm	80mm Ring
800040	Platinum 2 Piece Flat Baseplate	Pre-Cut 40mm	80mm Ring
800045	Platinum 2 Piece Flat Baseplate	Pre-Cut 45mm	80mm Ring
800050	Platinum 2 Piece Flat Baseplate	Pre-Cut 50mm	80mm Ring



Platinum 2PC Soft Convex 80mm Ring Size Baseplates

Units per box: 5

Product Code	Product Description	Size	Type
800015C	Platinum 2 Piece Soft Convex Baseplate	CTF 15-60mm	80mm Ring



Platinum 2PC 80mm Ring Size Closed Bags

Units per box: 30

Product Code	Product Description	Size
807415	Platinum 2 Piece Closed Bag	Standard
807615	Platinum 2 Piece Closed Bag	Maxi



Platinum 2PC 80mm Ring Size Drainable Bags

Units per box: 30

Product Code	Product Description	Size
818415	Platinum 2 Piece Drainable Bag	Standard
818615	Platinum 2 Piece Drainable Bag	Maxi



This is Maryrose... She's Family.

Maryrose is from Derry in Northern Ireland, and is the youngest of 10 children! Born with one kidney and an extra heart valve, she has been using Respond for 2 years.

"The service is great, reliable and always on time."

Joining Respond has never been so easy, get in touch today and let us do the rest!
#JoinTheRespondFamily



Phone
0800 220 300
0800 028 6848 (N.Ireland)



Online
respond.co.uk



Email
hello@respond.co.uk



Scan me
with your smartphone

Platinum with Vitamin E



Closed Bag.....	20
Soft Convex Closed Bag	20
Firm Convex Closed Bag	21
Drainable Bag.....	21
Soft Convex Drainable Bag.....	21
Firm Convex Drainable Bag	22
Soft Convex Urostomy Bag.....	22
Firm Convex Urostomy Bag	22
Soft Convex High Output Bag.....	23
Contour Soft Convex High Output Bag	23
Firm Convex High Output Bag	23
Contour Firm Convex High Output Bag	23
Drainable Post Op Bags	24
Platinum with Vit E Elements.....	24

Features & benefits

Infused with 1% Vitamin E

Promotes healthy skin by moisturising the delicate skin around the stoma. Starts to work in 15 minutes of application, and continues for at least 72 hours*.

Foam-backed baseplate

For added comfort.

Split fabric backing

For ease of use and filter access.

Available in flat, soft and firm convex.

High Output

High Output option available, with both contour and non-contour baseplates.

Urostomy Soft Bung outlet

Providing both security and comfort.

Discreet Velcro® outlet

Easy opening and positioning.



* In-vitro analysis undertaken using the High-Performance Liquid Chromatography (HPLC) technique

Platinum with Vitamin E Closed Bags

Units per box: 30

Sizes	Mini	Standard	Maxi
CTF 12-65mm	207212E	-	-
CTF 12-40mm	-	207412E	207612E
Pre-Cut 25mm	-	207425E	-
Pre-Cut 30mm	-	207430E	-
Pre-Cut 35mm	-	207435E	-
Pre-Cut 40mm	-	207440E	-
Pre-Cut 45mm	-	207445E	-
Pre-Cut 50mm	-	207450E	-
CTF 15-80mm	-	208415E	-



Platinum with Vitamin E Soft Convex Closed Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-65mm	407212E	-	-
CTF 12-40mm	-	407412E	407612E
Pre-Cut 25mm	-	407425E	-
Pre-Cut 30mm	-	407430E	-
Pre-Cut 35mm	-	407435E	-
Pre -Cut 40mm	-	407440E	-
CTF 15-60mm	-	408415E	-



Platinum with Vitamin E Firm Convex Closed Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-65mm	407212E-FC	-	-
CTF 12-40mm	-	407412E-FC	407612E-FC
CTF 15-60mm	-	408415E-FC	-



Platinum with Vitamin E Drainable Bags

Units per box: 30

Sizes	Mini	Standard	Maxi
CTF 12-65mm	218212E	218412E	218612E
Pre-Cut 25mm	-	218425E	-
Pre-Cut 30mm	-	218430E	-
Pre-Cut 35mm	-	218435E	-
Pre-Cut 40mm	-	218440E	-
CTF 15-80mm	-	220415E	-



Platinum with Vit E Soft Convex Drainable Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	418212E	418412E	418612E
Pre-Cut 25mm	-	418425E	-
Pre-Cut 30mm	-	418430E	-
Pre-Cut 35mm	-	418435E	-
Pre-Cut 40mm	-	418440E	-
CTF 15-60mm	-	420415E	-



Platinum with Vit E Firm Convex Drainable Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	418212E-FC	418412E-FC	418612E-FC
CTF 15-60mm	-	420415E-FC	-



Platinum with Vit E Soft Convex Urostomy Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	421212E	421412E	421612E
Pre-Cut 20mm	-	421420E	-
Pre-Cut 25mm	-	421425E	-
Pre-Cut 30mm	-	421430E	-
CTF 15-60mm		422415E	



Platinum with Vit E Firm Convex Urostomy Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	421212E-FC	421412E-FC	421612E-FC
CTF 15-60mm	-	422415E-FC	-



Platinum with Vit E Soft Convex High Output Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	-	-	318712E
CTF 15-60mm	-	-	320715E



Platinum with Vit E Contour Soft Convex High Output Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	-	-	418712E
CTF 15-60mm	-	-	420715E



Platinum with Vit E Firm Convex High Output Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
CTF 12-40mm	-	-	318712E-FC
CTF 15-60mm	-	-	320715E-FC



Platinum with Vit E Contour Firm Convex High Output Bags

Units per box: 30

Sizes	Mini	Standard	Maxi
CTF 12-40mm	-	-	418712E-FC
CTF 15-60mm	-	-	420715E-FC




Platinum with Vitamin E Drainable Post Op Bags

Sizes	Standard	Maxi	Units per box
CTF 12-65mm	218412E-PO	218612E-PO	30
CTF 15-40mm	418412E-PO	418612E-PO	10



Platinum with Vitamin E Elements

Product	Units per box	Product Code
Barrier Cream	1 tube (60g)	130100E
Barrier Film Spray	1 can (50ml)	130200E
Barrier Film Wipes	30 sachets	130300E
Adhesive Remover Spray	1 can (50ml)	130400E
Adhesive Remover Wipes	30 sachets	130500E



Platinum



Contour Closed Bag	27
Closed Bag.....	28
Contour Soft Convex Closed Bag.....	28
Soft Convex Closed Bag	29
Contour Drainable Bag	29
Drainable Bag.....	30
Contour Soft Convex Drainable Bag.....	30
Soft Convex Drainable Bag	31
Contour Soft Convex Urostomy Bag	31
Soft Convex Urostomy Bag.....	32
Drainable Post Op Bag.....	32

Features & benefits

Advanced filter system

The high-performance water-proof filter provides controlled airflow, which minimises the risks of ballooning and pancaking.

Soft Convexity

The Convex adhesive baseplate is soft, comfortable and flexible, and reduces leakages by conforming to the body.

Contour and non-contour

Platinum baseplates are available in contour and non-contour.

Soft and breathable fabric

Skin is able to breathe through the strong and durable non-woven fabric.

Discreet Velcro® outlet

Easy opening and positioning.

Foam-backed adhesive baseplate

Providing greater comfort and flexibility whilst maintaining superb adhesion.

Bung and Tap outlets

Platinum Urostomy bags are available with bung or tap option.



Platinum Contour Closed Bags

Units per box: 30

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-65mm	207212	207412	207612
Pre-Cut 25mm	-	207425	-
Pre-Cut 30mm	-	207430	-
Pre-Cut 32.5mm	-	207432	-
Pre-Cut 35mm		207435	
Pre-Cut 37.5mm		207437	
Pre-Cut 40mm		207440	
Pre-Cut 45mm		207445	
Pre-Cut 50mm		207450	
CTF 15-80mm		208415	



Pelican Ostomate Blogs

Our ambassadors regularly write blog posts with the aim of raising awareness and providing an insight into stoma life.



To read, scan the QR code or visit www.pelicanhealthcare.co.uk/blog

Platinum Closed Bags

Units per box: 30

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 12-65mm	107212	107112	107412	107312	107612	107512
Pre-Cut 25mm	-	-	107425	-	-	-
Pre-Cut 30mm	-	-	107430	-	-	-
Pre-Cut 35mm	-	-	107435	107335	-	-
Pre-Cut 40mm	-	-	107440	107340	-	-
Pre-Cut 45mm	-	-	107445	-	-	-
Pre-Cut 55mm	-	-	107455	-	-	-
CTF 15-80mm	-	-	108415	108315	-	-

Platinum Contour Soft Convex Closed Bags

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	407212	407412	407612
Pre-Cut 20mm	-	407420	-
Pre-Cut 25mm	-	407425	-
Pre-Cut 30mm	-	407430	-
Pre-Cut 35mm	-	407435	-
Pre-Cut 40mm	-	407440	-
CTF 15-60mm	-	408415	-



Platinum Soft Convex Closed Bags

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	307212	307412	307612
Pre-Cut 25mm	-	307425	-
Pre-Cut 30mm	-	307430	-
Pre-Cut 35mm	-	307435	-
CTF 15-60mm	-	308415	-



Platinum Contour Drainable Bags

Units per box: 30

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-65mm	218212	218412	218612
Pre-Cut 25mm	-	218425	-
Pre-Cut 30mm	-	218430	-
Pre-Cut 32.5mm	-	218432	-
Pre-Cut 35mm	-	218435	-
Pre-Cut 40mm	-	218440	-
Pre-Cut 45mm	-	218445	-
CTF 15-80mm	-	220415	-



Platinum Drainable Bags

Units per box: 30


Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-65mm	118212	118412	118612
Pre-Cut 25mm	-	118425	-
Pre-Cut 30mm	-	118430	-
Pre-Cut 35mm	-	118435	-
Pre-Cut 37.5mm	-	118437	-
Pre-Cut 40mm	-	118440	-
CTF 15-80mm	-	120415	-



Platinum Contour Soft Convex Drainable Bags

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	418212	418412	418612
Pre-Cut 20mm	-	418420	-
Pre-Cut 25mm	-	418425	-
Pre-Cut 30mm	-	418430	-
Pre-Cut 35mm	-	418435	-
Pre-Cut 40mm	-	418440	-
CTF 15-60mm	-	420415	-



Platinum Soft Convex Drainable Bags

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	318212	318412	318612
Pre-Cut 25mm	-	318425	-
Pre-Cut 30mm	-	318430	-
Pre-Cut 35mm	-	318435	-
CTF 15-60mm	-	320415	-



Platinum Contour Soft Convex Urostomy Bags with Bung Outlet

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	421212	421412	421612
Pre-Cut 25mm	-	421425	-
Pre-Cut 30mm	-	421430	-
CTF 15-60mm	-	422415	-



Platinum Contour Soft Convex Urostomy Bags with Tap Outlet

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	-	423412	423612
Pre-Cut 25mm	-	423425	-
Pre-Cut 30mm	-	423430	-
CTF 15-60mm	-	424415	-



Platinum Soft Convex Urostomy Bags with Bung Outlet

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	321212	321412	321612
CTF 15-60mm	-	322415	-



Platinum Soft Convex Urostomy Bags with Tap Outlet

Units per box: 10

Sizes	Mini Clear	Standard Clear	Maxi Clear
CTF 12-40mm	323212	323412	-
CTF 15-60mm	-	324415	-



Platinum Drainable Post Op Bags

Sizes	Mini	Standard Clear	Maxi Clear	Units per box
CTF 12-65mm	-	118412-PO	118612-PO	30
CTF 12-40mm	-	318412-PO	318612-PO	10

Select



Afresh Closed Bag.....	35
Closed Bag.....	36
Soft Convex Closed Bag	36
Afresh Drainable Bag	37
Soft Convex Drainable Bag	38
Urostomy Bag	40
Convex Urostomy Bag	41
Select Minuet Stoma Caps.....	42

Features & benefits

Coloured guides

Coloured measuring guide and thumb tab.

Soft Convexity

Available in soft convex.

Window

Clear and Opaque medical film options.

Split fabric backing

Fabric

Silky soft breathable fabric.

Urostomy


Extremely soft and flexible bung.

Size

Available in a wide range of sizes.

Drainable outlet available

Discreet and secure Velcro® outlet.



Select Afresh Closed Bags

Units per box: 30

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 20-65mm	102300	102320	102220	102120	-	-
CTF 20-65mm	-	-	-	-	102420	102400
Pre-Cut 22.5mm	-	-	102222	-	-	-
Pre-Cut 25mm	-	102325	102225	102125	-	-
Pre-Cut 27.5mm	-	-	102227	102127	-	-
Pre-Cut 30mm	-	-	102230	102130	-	-
Pre-Cut 32.5mm	-	102332	102232	102132	-	-
Pre-Cut 35mm	-	102335	102135	102235	-	-
Pre-Cut 37.5mm	-	-	102237	-	-	-
Pre-Cut 40mm	-	-	102140	102240	-	-
Pre-Cut 45mm	-	-	-	102245	-	-
Pre-Cut 50mm	-	-	-	102250	-	-

Select Closed Bags

Units per box: 30

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 20-65mm	100620	100600	-	-	-	-
CTF 20-65mm	-	-	100720	100700	-	-

Select Soft Convex Closed Bags

Units per box: 10

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 12-40mm	-	103112	103412	103312	103612	103512
CTF 15-60mm	-	-	104415	104315	-	-
Pre-Cut 20mm	-	-	-	103320	-	-
Pre-Cut 25mm	-	-	103425	103325	-	-
Pre-Cut 27.5mm	-	-	103427	103327	-	-
Pre-Cut 30mm	-	-	103330	103430	-	103530
Pre-Cut 32.5mm	-	-	103432	103332	-	-
Pre-Cut 35mm	-	-	103335	103435	-	103535
Pre-Cut 37.5mm	-	-	-	103437	-	-
Pre-Cut 40mm	-	-	103440	103340	-	103540

Select Afresh Drainable Bags

Units per box: 30

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 15-65mm	112415	112315	112215	112115	112615	112515
Pre-Cut 20mm	-	112320	112220	-	-	-
Pre-Cut 25mm	-	-	112225	112125	112625	112525
Pre-Cut 27.5mm	-	112327	112227	112127	-	112527
Pre-Cut 30mm	-	112330	112230	112130	112630	112530
Pre-Cut 32.5mm	-	-	112232	112132	-	-
Pre-Cut 35mm	-	-	112235	112135	-	112535
Pre-Cut 37.5mm	-	-	-	112137	-	-
Pre-Cut 40mm	-	-	112240	112140	112640	-
Pre-Cut 45mm	-	-	112245	112145	-	-
Pre-Cut 50mm	-	-	112250	112150	-	-

Select Soft Convex Drainable Bags

Units per box: 10

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 12-40mm	113212	113112	113412	113312	113612	113512
CTF 15-60mm	-	-	117415	117315	-	-
Pre-Cut 20mm	-	-	113420	113320	-	113520
Pre-Cut 22.5mm	-	-	113422	113322	-	113522
Pre-Cut 25mm	-	-	113425	113325	-	113525
Pre-Cut 27.5mm	-	-	113427	113327	-	113527
Pre-Cut 30mm	-	-	113430	113330	-	113530
Pre-Cut 32.5mm	-	-	113432	113332	-	113532
Pre-Cut 35mm	-	-	113435	113335	-	113535
Pre-Cut 37.5mm	-	-	113437	113337	-	-
Pre-Cut 40mm	-	-	113440	113340	-	113540



eakin
Cohesive®

Small but powerful change to a stoma care routine

eakin Cohesive® seals are a small change to your stoma care routine that can have a powerful impact every day.

*"I was probably about a year into my journey when I was introduced to the eakin Cohesive® seals and they've been amazing. They stretch, they don't lose how sticky they are, they fill gaps. For me it's not just a little bit of a change it has been a life change."
— Robyn, eakin Cohesive® user*


Scan me
to request
a sample



Select Urostomy Bags

Units per box: 30

Sizes	Mini		Standard		Maxi	
	Clear	Opaque	Clear	Opaque	Clear	Opaque
CTF 15-65mm	114415	-	114215	-	114615	-
Pre-Cut 20mm	-	-	114220	-	-	-
Pre-Cut 22.5mm	-	-	114225	-	-	-
Pre-Cut 27.5mm	-	-	114230	-	-	-
Pre-Cut 30mm	-	-	114232	-	-	-
Pre-Cut 32.5mm	-	-	114235	-	-	-
Pre-Cut 40mm	-	-	114245	-	-	-



Select Soft Convex Urostomy Bags

Units per box: 10

Sizes	Mini	Standard	Maxi
	Clear	Clear	Clear
CTF 12-40mm	115212	115412	115612
CTF 15-60mm	-	119415	-
Pre-Cut 20mm	-	115420	-
Pre-Cut 22.5mm	-	115422	-
Pre-Cut 25mm	-	115425	-
Pre-Cut 27.5mm	-	115427	-
Pre-Cut 30mm	-	115430	-
Pre-Cut 32.5mm	-	115432	-
Pre-Cut 35mm	-	115435	-
Pre-Cut 40mm	-	115440	-

Select Minuet Stoma Caps

Units per box: 30

Product	Product Code
Select Minuet Stoma Cap (10-40mm) Clear	102380
Select Minuet Plus Stoma Cap (10-55mm) Clear	102382

Features and benefits

- Unique, comfortable and secure skin protector
- Silky soft fabric for comfort and discretion
- Split fabric backing, to allow access to cover the filter when required
- Hydrophobic filter system
- Available in two sizes – Select Minuet & Minuet Plus
- Ideal for irrigation, sports and leisure, and intimate moments



Neonatal and Paediatric



Select Afresh Paediatric	45
Select Afresh Neonatal	45
Select Soft Convex Closed Paediatric	45
Select Soft Convex Drainable Paediatric	46
Select Urostomy Paediatric	46
Select Paediatric.....	47
Classic Paediatric/Neonatal	47

Features & benefits

Waterproof filter

Designed to help reduce pancaking and ballooning.



Clipless system

Quick and easy to clean, with a secure Velcro® tab.

Ultra-thin, hypoallergenic skin protector

Moulds to the fit of the body, allowing for extra security.



Available in Urostomy

Bung gives ability to control flow when emptying.

Gentle against the skin

Soft fabric with Teddy Bear print.

Select Afresh Paediatric Bag with Filter (Velcro® closure)

Units per box: 30

Size	Fabric	Film	Product Code
CTF 10-50mm	Teddy Print	Clear	112701
CTF 10-50mm	Plain	Opaque	112702



Select Afresh Neonatal Bag with Filter (Velcro® closure)

Units per box: 30

Size	Fabric	Film	Product Code
CTF 0-40mm	Teddy Print	Clear	112703
CTF 0-40mm	Plain	Opaque	112704



Select Soft Convex Closed Paediatric Bag with Filter (Velcro® closure)

Units per box: 10

Size	Fabric	Film	Product Code
CTF 5-30mm	Teddy Print	Clear	112709
CTF 5-30mm	Plain	Clear	112711



Select Soft Convex Drainable Paediatric Bag with Filter (Velcro® closure)

Units per box: 10

Size	Fabric	Film	Product Code
CTF 5-30mm	Teddy Print	Clear	112708
CTF 5-30mm	Plain	Clear	112710



Select Urostomy Paediatric Bags

Size	Fabric	Film	Units per box	Product Code
CTF 10-50mm	Teddy Print	Clear	30	112706



Select Convex Urostomy Paediatric Bags

Size	Fabric	Film	Units per box	Product Code
CTF 5-30mm	Teddy Print	Clear	10	112712
CTF 5-30mm	Plain	Clear	10	112713



Select Paediatric Bag with Filter (Clip closure)

Units per box: 30

Size	Fabric	Film	Product Code
CTF 10-50mm	Teddy Print	Clear	101608

Select Paediatric Bag (Clip closure)


Units per box: 30

Size	Fabric	Film	Product Code
CTF 10-50mm	Plain	Clear	101600
CTF 10-50mm	Plain	Opaque	101602

Classic Paediatric/Neonatal Bag (Clip closure)

Units per box: 30

Size	Description	Film	Product Code
CTF 7-40mm	Classic Drainable	Clear	120705



Meet our Pelican Ambassadors

We are so grateful to work with our amazing group of ambassadors. Our inspirational ostomates provide insight, support and raise awareness through discussion groups, idea sessions, blog posts and by creating engaging content. Keep an eye out for these familiar faces!



JACK
@ostomyjack



LOUISE
@crohnsfighting



ANGELINA
@angelinatruman



ANT
@ibdlife



KEITH
@keiththom2014



CHARLOTTE
@char.crohnsftwilson



LAURA
@2bagsforlife



RACHEL
@gutsy.mum



SUMMER
@summers_stoma

Why not give them and us a follow!



@pelicanhealthcare



PelicanHealthcareLtd



@pelicanhealth



Pelican Healthcare Ltd

Elements



Platinum with Vit E Elements...	51
Select Minuet Stoma Caps.....	51
2-litre Night Drainage Bag	51
Pelican Elements	52
Support Belts.....	53

Features & benefits

The Pelican Elements range provides a wide selection of accessory products designed to solve a range of problems that you and your patients may encounter.

The range includes the Platinum with Vitamin E Elements, infused with 1% Vitamin E, designed to not only look after your patients general stoma care but also the healthy skin surrounding it.



Select Minuet Stoma Caps

Product	Units per box	Product Code
Select Minuet Stoma Cap (10-40mm) Clear	30 per pack	102380
Select Minuet Plus Stoma Cap (10-55mm) Clear	30 per pack	102382



Pelican 2-litre Night Drainage Bag

Product	Units per box	Product Code
Non-Drainable	10	133100
Drainable	10	133200



Platinum with Vitamin E Elements

Product	Units per box	Product Code
Barrier Cream	1 tube (60g)	130100E
Barrier Film Spray	1 can (50ml)	130200E
Barrier Film Wipes	30 sachets	130300E
Adhesive Remover Spray	1 can (50ml)	130400E
Adhesive Remover Wipes	30 sachets	130500E

Pelican Elements

Product	Units per box	Product Code
Pelican Protect Plus non-sting barrier film wipes	30 sachets	130260
Pelican Protect Plus non-sting barrier film spray	1 can (50ml)	130261
Pelican Contour Flange Extenders	30	130860
Pelican Perform sachets	60 sachets	130160
Pelican Release non-sting adhesive remover wipes	30 sachets	130360
Pelican Release non-sting adhesive remover spray	1 can (50ml)	130460
Pelican Dor pouch deodorant	1 bottle (7 ml)	130103
Adjustable ostomy belt	1	133006




eakin® Support Belts 20cm Width

Description	Measurements	Product Code
Small	75-90cm	839051
Medium	90-105cm	839052
Large	105-120cm	839053
Extra Large	120-135cm	839054



eakin® Support Belts 26cm Width

Description	Measurements	Product Code
Small	75-90cm	839061
Medium	90-105cm	839062
Large	105-120cm	839063
Extra Large	120-135cm	839064



Are you looking for a reliable and supportive prescription service?



Hear what being part of the Respond service means to our ambassadors by scanning the QR code or visiting www.respond.co.uk/join-the-family/

Scan me with
your smartphone to
watch the video

It's a Family Thing.

Joining Respond has never been so easy, get in touch today and let us do the rest!
#JoinTheRespondFamily



Phone

0800 220 300

0800 028 6848 (N.Ireland)



Online

respond.co.uk



Email

hello@respond.co.uk



Providing NHS services

eakin® Seals



eakin freeseal ®	57
eakin Cohesive ® seals	57
eakin Cohesive ® paste	57
eakin Cohesive ® skin barriers ...	57

Features & benefits

eakin Cohesive® large

Gives the greatest coverage to smooth uneven skin surfaces such as dips or creases.

eakin Cohesive® slims

The leak preventer, for ongoing skin protection and confidence for stomas under 50mm.

eakin Cohesive® small

The extra protector, eakin®'s thickest seal for stomas under 50mm with frequent loose output.

eakin Cohesive® stomawrap

The time saver. Oval-shaped for larger or irregular shaped stomas over 50mm and post operative swelling. Unique wrap around - easy to use for limited dexterity or visual impairment.

eakin freeseal®

eakin®'s thinnest seal, at 1.8mm thin, ideal for convex drainable bags.



eakin freeseal®

Product Name	Size	Units per box	Product Code
eakin freeseal®	Diameter - 48mm Depth - 1.8mm	10	839007



eakin Cohesive® seals

Product Name	Size	Units per box	Product Code
eakin Cohesive® large	Diameter - 98mm Depth - 2.1mm	10	839001
eakin Cohesive® small	Diameter - 48mm Depth - 4.2mm	30	839002
eakin Cohesive® slims	Diameter - 48mm Depth - 2.8mm	30	839005
eakin Cohesive® stomawrap	Diameter - 85mm Depth - 3.0mm	10	839006



eakin Cohesive® paste

Product Name	Units per box	Product Code
eakin Cohesive® paste	60 grams	839010



eakin Cohesive® skin barriers

Product Name	Size	Units per box	Product Code
eakin Cohesive® skin barrier large	10 x 20cm	5	839003
eakin Cohesive® skin barrier small	10 x 10cm	5	839004





#BeThe
Change



Insights
Community

Would you like to play a part in helping us develop our next generation of products and services?

Our vision is to support ostomates in leading a full and dignified life – but we need your help to do this!

Our new insights community is all about giving you a voice and a platform to tell us what is important to you and your patients.

You can choose to be involved in as little or as much as you would like – it's completely up to you!



Insights - Community - Innovation

Please register your interest by scanning the QR code.



#BeTheChange



Eakin
Healthcare

eakin® Wound

eakin Wound Pouches™ with Bung	61
eakin Wound Pouches™ with Fold & Tuck....	62
eakin Access Windows™	62
eakin Cohesive ® skin barriers	62



Features & benefits

- A range of shapes and sizes to ensure suitability
- Enclosed tracing template for maximum accuracy
- Two closure methods (Bung and Fold & Tuck)
- Transparent film for viewing without pouch removal
- Ability to add an eakin Access Window™ to provide instant and repeated access to the wound
- Can be used in conjunction with eakin **Cohesive**® skin barriers, paste and seals to protect the skin
- Can be connected to continuous drainage



eakin Wound Pouches™ Bung Closure

Product Name	Size	Units per box	Product Code
Small	45 x 30mm	10	839260
Small Plus	86 x 60mm	10	839267
Medium	110 x 75mm	10	839261
Large	175 x 110mm	10	839262
Large (for vertical incision wounds)	250 x 80mm	5	839270
Extra Large (for horizontal wounds)	245 x 160mm	5	839263
Extra Large with Access Window (for horizontal wounds)	245 x 160mm	5	839264
Extra Large with Access Window (for vertical wounds)	245 x 160mm	5	839268
Extra Large (for vertical wounds)	245 x 160mm	5	839265
Extra Large (for vertical incision wounds)	290 x 130mm	5	839266



eakin Wound Pouches™ Fold & Tuck Closure

Product Name	Size	Units per box	Product Code
Small	45 x 30mm	10	839250
Medium	110 x 75mm	10	839251
Large	175 x 110mm	10	839252
Extra Large	245 x 160mm	5	839253



eakin Access Windows™

Product Name	Size	Units per box	Product Code
Large Access Window	125 x 75mm	5	839280
Small Access Window	60 x 60mm	5	839281



eakin Cohesive® skin barriers

Product Name	Size	Units per box	Product Code
eakin Cohesive ® skin barrier large	10 x 20cm	5	839003
eakin Cohesive ® skin barrier small	10 x 10cm	5	839004



To order stock or for any additional information:

Call **0800 318 282**

Visit **pelicanhealthcare.co.uk**

Email **contactus@pelicanhealthcare.co.uk**

 PelicanHealthcare

 PelicanHealthcareLtd

 PelicanHealth

 Pelican Healthcare Ltd

PELICANPORTFOLIO/2



© 2023 Pelican Healthcare Ltd

pelican
HEALTHCARE



Eakin
Healthcare

Working together. Improving lives.